


DIXION OPEN REANIMATION SYSTEMS

DIXION Open Reanimation Systems with three control modes (pre-warm mode, manual mode and baby mode) are controlled by a micro-computer. Set temperature and baby temperature can be displayed separately.

- X-ray cassette under the bassinet
- APGAR timer
- Phototherapy time timer
- Failure alarm indication



	Baby Guard W-1145	Baby Guard W-1140
Servo control t° air	•	•
Servo control t° skin	•	•
Heater "Radiant heat"	•	•
Adjustable bassinet inclination	•	•
Adjustable angle of phototherapy/warmer module	•	•
Reclining side panels	•	•
X-ray cassette holder	•	•
LED phototherapy	•	•
LCD display with backlight	•	•
Data storage	•	•
Lower bed heating	•	•
Gel mattress	•	○
Stand height adjustment	•	○
Built-in battery	•	•
2 backlight modes in- use and night	•	•
Case for accessories	•	•
Apgar timer	•	•

○ — optional